


Quality through Experience

ERWEKA®

Certificate

Interchangeability of ERWEKA vessel

ERWEKA applies USP compliant round-bottomed, numbered 1000ml vessels as part of the dissolution assembly. We herewith confirm that our vessels numbered from Q1001 are interchangeable and are versatile usable in every test position. Measurements with regard to vessel verticality, centricity and height of stirring elements haven proven full vessel compatibility and compliance with the FDA requirements "Mechanical Calibration".



Heusenstamm, 2011-03-15

ppa Eggert

ppa. Dr. Renate Eggert
ERWEKA Head of Quality

ERWEKA GmbH
Ottostr. 20-22
D-63150 Heusenstamm

Telefon: (06104) 69 03 0
Telefax: (06104) 69 03 40

Geschäftsführer
Managing Director
Werner G. Müller
Christian G. Müller
Amtsgericht Offenbach
HRB 2382