

Konformitätserklärung

Certificate of Compliance

Prüfkorb für Zerfallzeittester nach USP/EP

Basket rack for Disintegration Tester according to USP/EP

Prüfkorb:

Basket rack:

Material Distanzbolzen: 1.4571^①

Material of distance rods:

Außendurchmesser: 97 mm

Outside diameter:

Röhrchen:

Tubes:

Material: Glas

Material: Glass

Innendurchmesser: 33,0 ± 0,5 mm


Inside diameter:

Wanddicke: 2,5 ± 0,5 mm

Thickness of the wall:

Länge: 77,5 ± 2,5 mm

Length:



^① German stainless steel according to EN 10027-2,
similar to AISI/SAE 316Ti

Sieb:**Sieve:**Material: 1.4401^②

Material:

Maschenweite: 2,0 ±0,2 mm

Mesh size:

Drahtdurchmesser: 0,63 ± 0,03 mm (entspricht Mesh 10)

Wire diameter:

(equivalent to mesh 10)

Schwimmdüse:

siehe separate Konformitätserklärung

Disc:

see separate Declaration of Conformity

Für den Prüfkorb unserer Zerfallzeittester bescheinigen wir in alleiniger Verantwortung die Übereinstimmung mit den Forderungen der 9thUSP41, Dietary Supplements Chapter <2040> und der EP 9th Edition (Kapitel 2.9.1, Apparatur B).

For the basket rack of our disintegration tester we declare in our sole responsibility the conformity with the 9thUSP41, Dietary Supplements Chapter <2040> and the EP 9th Edition (Chapter 2.9.1, apparatus B).

Langen, 2019-05-20

i.A.



Sevgi Klier

- Head of Integrated Management System -

^② German stainless steel according to EN 10027-2, similar to AISI/SAE 316